



RNase A Contamination: Unpredictable and Catastrophic

*Prevent Contamination in your Lab with High
Irradiance UV LED*

A Phoseon Technology White Paper
January 2019

Ribonuclease (RNase) Contamination: Unpredictable and Catastrophic

RNases are an ongoing problem for experiments requiring full length RNA. Even trace amounts of RNase have a big impact on RNA sequencing.

The single most important aspect of RNA protocols is isolating and maintaining full length, un-degraded RNA for analysis or use as a reaction substrate. Whether preparing total RNA libraries for NGS or looking at individual RNAs (iCLIP), un-degraded RNA is pivotal. Hindering these processes is RNase.

Degradation by nucleases such as RNase is a recurring laboratory handling issue which requires diverse cleaning methods that are both time-consuming and costly. Repeated cleaning not only leaves residues, another source of contamination, but can also affect the integrity of equipment (corroding metal and degrading plastics). It is also not always clear when these cleaning methods should be implemented.

Furthermore, repeated exposure to cleaning solutions or soaking may corrode metal or degrade plastic surfaces. How clean is clean enough? Clean enough occurs when you don't need to repeat lengthy protocols because of degraded RNA. Even trace amounts of RNase have a big impact on RNA sequencing, due to its catalytic action. Not clean leads to time and money loss.

This white paper makes the case for why UV LED technology deserves serious consideration by RNA sequencing labs for controlling ribonuclease. It describes Phoseon's findings related to LED light engines for the inactivation of RNases in a laboratory setting.



Quantification of Ambient RNase A Contamination

RNA Institutional knowledge says the level of contamination is directly proportional to time. Thus one could wait a specified amount of time to clean surfaces or equipment. However, this approach does not address the seemingly random contamination that often plagues RNA-work, a particular concern when it comes to microplates.

The more time microplates spend removed from their sterile environment the greater the contamination. Unsurprisingly, we found that a microplate exposed to a business cubicle environment for 1.5 months showed higher levels of contamination overall than one exposed for just five minutes (figure 1A-B). But contamination is not only proportional to total time of exposure. It also occurs sporadically, resulting in unpredictable variations in contamination levels between individual wells of a microplate.

Contamination occurs sporadically, resulting in unpredictable variations in contamination levels between individual wells of a microplate.





Figure 1: Long-term vs short-term ambient contamination exposure
A flat-bottomed, black opaque microplate was exposed to an office cubicle environment for 1.5 months (a) and a second microplate was exposed to an office lunch-room environment for 5 minutes (b). Sample wells were suspended and fluorometrically assayed for RNase activity (RNase Alert IDT, Gemini XPS).

In our studies of RNase A inactivation with a high intensity UV LED scanning array, we found the burden RNase contamination can be significantly reduced.

The microplate exposed to a lunchroom environment for five minutes showed one well with significant levels of contamination above the negative control while all other wells showed no significant contamination (figure 1B). This one well represents instantaneous contamination irrespective of the total time of exposure. Implying that simply limiting exposure to environmental RNase does not ensure your microplate is free of contamination. Up to this point laboratories doing RNA work have been forced to bear the weight of this random ambient contamination in their microplates, repeating lengthy protocols while hoping the contamination is not also repeated. In our studies of RNase A inactivation with a high intensity UV LED scanning array, we found the burden RNase contamination can be significantly reduced.

Working from previous knowledge that 275nm and 365nm sources work synergistically to speed RNase A inactivation (1), we used these wavelengths on Greiner Bio-One flat-bottomed black opaque microplates. Sample wells were first treated with 0.01u/mL of RNase A. We then used 1.4 W/cm² of 275nm light and 700 mW/cm² of 365nm light (figure 2C) on a scanning array to decontaminate the doped wells. Once efficacy was confirmed on Greiner plates, we began testing microplates of other materials, colors, and well shapes.

Each plate was first melt-tested to determine the intensity of UV light it could withstand. Some plate materials proved to be incompatible with the high-intensity UV at RNase inactivation wavelengths (figure 2B).



Why UV LED Technology?

*Inactivate RNase A with High-intensity
UV LED Light Engines*



We then determined the minimum intensity of light necessary to successfully inactivate RNase A in each plate type. Figure 2 outlines the various kinds of plates that were tested and what levels of UV were necessary.

With a high-intensity, scanning UV LED array we effectively decontaminated the most common microplate types in less than five minutes. This protocol does not leave any residue as chemicals and sprays often do. Furthermore, by implementing this protocol prior to every run, the integrity of your RNA library or sequence is secure.

With a high-intensity, scanning UV LED array we effectively decontaminated the most common microplate types in less than five minutes.



KeyPro™ KP100 Decontamination System

High-intensity LED Light Engines Effectively Inactivate RNase A

High-intensity UV light has been shown to rapidly, effectively, irreversibly inactivate RNase.

Working with RNA can be intimidating. Environmental RNase contamination sources include microbial contamination from room air as well as RNases from human skin, hair, or saliva. RNase inactivation methods range from DEPC treatment followed by autoclaving to more involved methods such as: chemical decontamination of surfaces, baking glassware, rinsing equipment in RNase-free water after chemical treatment, and frequent glove changes - all while continually using freshly opened disposables. Such cleaning methods can be costly in terms of money but more importantly they are time consuming, slowing research throughput and likely leading to erroneous results. Now there is a better solution.

Phoseon Technology is the first to develop a UV LED system that surpasses 2.5 W/cm², significantly higher than the levels reached by other technologies in the market. This high-intensity UV light has been shown to rapidly, effectively, irreversibly inactivate RNase. This milestone development provides scientists, researchers and equipment manufacturers the capability to rapidly and reliably control RNase contamination.

RNA Protocols that benefit from UV inactivation of RNase:

- 1) *Ultra-low input and Single-cell RNA sequencing*
- 2) *Ribosome profiling*
- 3) *RNA Exome Capture sequencing*
- 4) *Targeted RNA sequencing*
- 5) *Small RNA sequencing*
- 6) *Total RNA sequencing*
- 7) *mRNA sequencing*
- 8) *CRAC (Crosslinking And cDNA analysis)*
- 9) *iCLIP (individual-nucleotide resolution Cross-Linking and ImmunoPrecipitation)*
- 10) *NGS of RNAs*

Why UV LED Technology?

*Inactivate RNase A with High-intensity
UV LED Light Engines*



Conclusion

- Specific wavelengths, and wavelength combinations, of UV LED irreversibly inactivate RNase A when used at a sufficiently high irradiance.
- The high irradiance necessary for RNase A inactivation is made possible by Phoseon's SLM technology.
- Wavelengths of 275 nm and 365 nm interact synergistically resulting in faster inactivation, at lower irradiances, than is achievable with either wavelength alone.
- UV LED inactivation of RNase A is much faster than conventional methods and does not leave any chemical residue on surfaces.

In short, RNases are an on-going problem for experiments requiring full length RNA. Application of UV LED technology can benefit researchers through improved reliability of starting materials, shorter time required for preparation and inactivation of RNases, all while protecting valuable RNA samples from degradation and chemical contamination.



Author

Theresa Thompson, PhD Application Scientist, Phoseon Technology

For more information about Phoseon Technology products and services, please contact:

Theresa Thompson, theresa.thompson@phoseon.com

References

1. Irreversible inactivation of RNase A on a surface by UV LED; Thompson T, Eliason G, Pasquantonio J; American Society of Human Genetics Scientific Sessions 2017 (poster)
2. Rapid inactivation of RNase A by high irradiance UV LEDs; Thompson T, Eliason G, Pasquantonio J; American Society of Cell Biology/European Molecular Biology Organization Scientific Sessions 2017 (poster)
3. The Importance of RNase Inactivation in RNA Sequencing Labs; Thompson T, Taggard K, Pasquantonio J; Whitepaper
4. Avoiding Ribonuclease Contamination Tools and Guidelines; Biolabs Website
5. Avoiding RNase contamination; Wikipedia
6. Is DEPC treatment of glassware and plasticware necessary for RNA extraction?; ResearchGate

Why UV LED Technology?
*Inactivate RNase A with High-intensity
UV LED Light Engines*



About Phoseon Technology

The world leader since 2002, Phoseon Technology pioneered the use of LED technology for Life Science and Industrial Curing applications. Phoseon delivers innovative, highly engineered, patented LED solutions. The company is focused 100% on LED technology and provides worldwide support.

Contacts

For more information about Phoseon Technology suite of products, visit <http://www.phoseon.com/> or call (503) 439-6446

Copyright © 2018 by Phoseon Technology. All rights reserved.

The information contained herein has been obtained by Phoseon Technology from sources believed to be reliable. However, because of the possibility of human or mechanical error by these sources, Phoseon Technology does not guarantee the accuracy, adequacy, or completeness of any information and is not responsible for any errors or omissions or for the results obtained from use of such information.

This material is the property of Phoseon Technology or is licensed to Phoseon Technology. This material may only

be reproduced, transmitted, excerpted, distributed, or commingled with other information, with the express written permission of Phoseon Technology. The user of this material shall not edit, modify, or alter any portion.

Requests for use may be submitted in writing to www.phoseon.com/contact. Any material quoted from this publication must be attributed to "UV Power or Irradiance? Understanding the Difference and Why Both Matter," published by Phoseon Technology, © 2018 by Phoseon Technology. All Rights Reserved."